




GUIDE TO YOUR NEW AVACTA ALLERGY+ ENVIRONMENTAL TEST RESULTS

Provision of the AU value enables the determination of more subtle differences in reactivity; this may help guide in the selection of allergens for inclusion in immunotherapy (ASIT).


Coloured bars offer clear visual comparison between allergen reactivity and are a great visual tool when discussing the results with owners.

Despite frequently showing elevated IgE levels, moulds are seldom clinically significant and results should be interpreted with this in mind.



Allergens highlighted in green have test values higher than the vertical cut-off line and should be considered for allergen-specific immunotherapy (ASIT) and allergen avoidance, providing the correct diagnostic work-up has been followed and relevant clinical signs are present.

Guidelines to aid interpretation. If you have any questions please contact our Technical Support Team. If it is more convenient, please use the 'Contact Us' form on the website where you can book an allocated call back slot at a time convenient for you. Remember, you can also access 24/7 support via our Practice Portal at avactaanimalhealth.com/login.



For further information or assistance, please contact our Technical Support Team on **0800 3 047 047** or at technical.support@avacta.com

Avacta Animal Health Ltd, Unit 651, Street 5, Thorp Arch Trading Estate, Wetherby, UK, LS23 7FZ
T – 0800 3 047 047 E – technical.support@avacta.com W – avactaanimalhealth.com